

Introduction

The capture and molecular analysis of cell-free DNA (cfDNA) released from cells under certain pathological conditions, including cancer, provides a non-invasive alternative to currently used practices via use of liquid biopsy samples.

The rapid development of next-generation sequencing (NGS) technologies in recent years, has led to a significant reduction in sequencing cost with improved accuracy.

In liquid biopsy research, NGS can be applied to sequence circulating tumour DNA (ctDNA; DNA fragments released by tumour cells). As such, NGS could provide molecular characterisation of tumour-specific genomic alterations in ctDNA, allowing for non-invasive and real-time monitoring of disease development.

The aim of this study is to evaluate a modified version of the OGT™ Universal NGS Workflow Solution in conjunction with a custom SureSeq™ targeted panel for use with cfDNA.

Methods

A modified version (not commercially available) of the OGT Universal NGS Workflow Solution was used throughout this study (Fig 1). The approach offers a streamlined NGS library preparation protocol with Unique Molecular Identifiers (UMIs) and unique dual indexing (UDI), followed by hybridisation-based target enrichment.

For this capability study, we used a 64 kb SureSeq myPanel™ Custom panel comprising 213 target exons in 40 genes.

Sequencing was conducted using either MiSeq® V2 300 or NextSeq® 500/550 Mid output chemistry (Illumina).




Figure 1: Modified OGT SureSeq workflow for cfDNA. The workflow allows users to go from extracted cfDNA to sequencer in 1.5 days with minimal handling time.

Samples

We tested a commercially available cfDNA standard (SeraSeq® ctDNA Complete™ Reference Material AF1%, SeraCare) with 10–50 ng DNA input. To mimic different variant allele frequencies (VAF), we diluted the reference DNA to create samples with a frequency range of 0.5 – 1%.

Bioinformatic Analysis

Sequencing data analysis was performed using OGT's proprietary Interpret NGS Analysis Software. Specifically, the software was used for read alignments, UMI deduplications, coverage calculations and variant calling. The data can also be visualised within the software using an IGV (Integrated Genomic Viewer).

Results I

Generation of high-quality sequencing libraries from low input cfDNA

We have successfully generated high-quality sequencing libraries using 10 – 50 ng cfDNA.

Hybridisation-based enrichment of these libraries has achieved a high unique coverage of ≥1,000x from as little as 10 ng cfDNA input (minimum 11 M reads).

We have achieved very high mean coverages (1,500x – 2,000x) from 25 ng and 50 ng DNA inputs (minimum 12 M reads).




Figure 2: Mean Target Coverage for libraries generated from 10, 25 or 50 ng cfDNA. Error bars reflect standard deviation of replicates (n=5).

Results II

Highly uniform coverage allowing confident SNVs detection in cfDNA samples

The hybridisation-based OGT Universal workflow coupled with OGT's expert bait design deliver excellent depth of coverage and coverage uniformity for 10 – 50 ng cfDNA (Fig 3), allowing consistent detection of SNVs (Fig 4).




Figure 3: Illustration of the excellent coverage uniformity of BRCA1 exon 10 in cfDNA samples with 10, 25 or 50 ng input amount. Depth of coverage per base (grey). Gene coding region as defined by RefSeq (bottom track)




Figure 4: Example of BRCA1 exon 10 missense variant Pro871Leu (rs799917) with frequency of 20% identified in 10 ng cfDNA. Shown are the genomic location (top), depth of coverage per base (grey bars), aligned example reads (grey lines) and the reference sequence.

Results III

For proof of principle experiments, we tested 0.5% VAF for inputs of 25 and 50 ng, and 1% VAF for 10 ng of SeraCare cfDNA reference material. These VAF were chosen to match a theoretical detection limit of ≥10 supporting reads for variants based on the Mean Target Coverage we achieved for each input in Results I.

Confident detection of SNVs and indels with 0.5% VAF using 25 ng or 50 ng cfDNA

Gene	Variant	25 ng 0.5% VAF		50 ng 0.5% VAF	
		Total reads	Variant reads	Total reads	Variant reads
AKT1	c.49G>A	3,998	48	7,892	67
APC	c.4348C>T	6,989	48	10,497	52
BRAF	c.1799T>A	4,765	33	7,096	60
CTNMB	c.121A>G	7,342	80	10,225	101
FGFR3	c.746C>G	4,643	51	8,710	71
KIT	c.2447A>T	4,087	28	5,618	57
KRAS	c.35G>A	4,809	32	6,613	82
NRAS	c.182A>G	5,033	42	6,800	40
PDGFRA	c.2525A>T	4,919	30	7,426	85
PIK3CA	c.3140A>G	4,940	22	7,147	35
RET	c.2753T>C	4,074	39	5,768	39
TP53	c.524G>A	7,234	69	10,790	105
TP53	c.743G>A	7,032	50	10,821	97
TP53	c.818G>A	4,821	40	6,498	59
PDGFRA	c.1694_1695insA	4,747	61	6,882	43
PIK3CA	c.3204_3205insA	3,522	19	5,222	28
APC	c.4666dup 6A>7A	3,042	7	4,240	17
TP53	c.723del	5,321	36	8,298	65
TP53	c.267del	3,246	27	4,707	57

A modified version of OGT's Universal workflow together with a SureSeq myPanel Custom panel successfully detected all anticipated SNVs (n=14) and Indels (n=5) with 0.5% VAF using 25 or 50 ng SeraCare cfDNA (Table 1).

The data also shows that with deeper sequencing, Mean Target Coverage of 5,400x and 7,700x can be achieved for 25 ng (29 M reads) and 50 ng (41 M reads) inputs, respectively. This suggests a VAF detection limit below 0.5% might be possible.

Table 1: Detection of SNVs and Indels in 25 and 50 ng SeraCare cfDNA reference standards with 0.5% VAF. Variant reads are filtered for unique UMIs only, and not for UMI family sizes.

Confident detection of SNVs and indels with 1% VAF using 10 ng cfDNA

Gene	Variant	10 ng 1% VAF	
		Total reads	Variant reads
AKT1	c.49G>A	3,909	34
APC	c.4348C>T	6,219	30
BRAF	c.1799T>A	4,543	37
CTNMB	c.121A>G	6,127	55
FGFR3	c.746C>G	2,404	21
KIT	c.2447A>T	3,601	27
KRAS	c.35G>A	3,654	46
NRAS	c.182A>G	4,455	36
PDGFRA	c.2525A>T	4,609	54
PIK3CA	c.3140A>G	4,382	45
RET	c.2753T>C	3,837	34
TP53	c.524G>A	5,196	74
TP53	c.743G>A	6,064	60
TP53	c.818G>A	4,249	67
PDGFRA	c.1694_1695insA	4,210	51
PIK3CA	c.3204_3205insA	3,330	27
APC	c.4666dup 6A>7A	2,986	30
TP53	c.723del	5,117	55
TP53	c.267del	2,085	16

We have achieved detection of all anticipated SNVs (n=14) and Indels (n=5) with 1% VAF using 10 ng SeraCare cfDNA (Fig 5, Table 2). A Mean Target Coverage of 4,200x was achieved with 43 M reads.




Figure 5: Three examples of detected SNVs with 1% VAF in 10 ng - PIK3CA exon 21 His1047Arg (COSM775), NRAS exon 3 Gln61Arg (COSM584) and KRAS exon 2 Gly12Asp (COSM521). Shown are the genomic location (top), depth of coverage per base (bars) and the reference sequence.

Table 2: Detection of SNVs and Indels in 10 ng SeraCare cfDNA reference standards with 1% VAF. Variant reads are filtered for unique UMIs only, and not for UMI family sizes.

Conclusions

- We have demonstrated that a modified version (not commercially available) of OGT Universal NGS Workflow Solution is suitable for very low inputs of cell-free DNA.
- High coverage uniformity after hybridisation-based enrichment using a custom panel allowed simultaneous detection of SNVs and indels in hundreds of target exons.
- We have shown that our approach reliably detected somatic mutations down to 0.5% VAF.
- These results confirm the capability of utilising this type of approach for cell-free DNA research applications.

