

# CytoSure



A Sysmex Group Company

## Disease-Focused Research Arrays



## Features

### **Targeted content and high-quality probe design**

- Detect single or multiple exonic CNVs in the genes that matter

### **Customisable multiplex designs**

- High-throughput, Cost-effective analysis

### **Simplified data interpretation**

- Using industry-leading CytoSure Interpret Software

### **Optimised labelling kits and integrated sample tracking probes**

- Confident analysis and reporting

### Reliable detection of copy number changes for a range of genetic disorders

CytoSure® Disease-Focused Research arrays are exon-focused microarray designs containing probe sets targeted to a specific disease, or groups of disease (Table 1). These disease-focused research arrays enable single-exon resolution copy number variation (CNV) analysis, allowing the detection of microdeletions and microduplications within targeted genes. The multiplex format enables a high-throughput approach to be taken, focused on genes with known relevance to the disorder — helping to minimise variants of unknown significance — perfectly complementing targeted next generation sequencing (NGS) panels or replacing multiple multiplex ligation-dependent probe amplification (MLPA) probe mixes.

### Targeted content

Inherited genetic disorders can be caused by a variety of chromosomal aberrations, including point mutations and small CNVs (Figure 1). Different methodologies are combined to accurately detect these changes, with one of the most successful combinations being targeted NGS and array comparative genomic hybridisation (aCGH)<sup>1</sup>. Using highly-targeted, exon-focused arrays as part of this combined approach has been shown to detect small CNVs (Figure 2).

The concept of using highly targeted gene-focused arrays to screen specific loci has been developed in collaboration with leading molecular genetics experts at Emory University and has now been implemented in multiple labs worldwide.




Figure 1: The mutation spectrum of molecular disorders includes point mutations and deletions and duplications. The prevalence of each type of variant is highly dependent on the disease studied.




Figure 2: CytoSure Interpret Software enables easy detection and visualisation of small aberrations. Shown here is a 801bp micro-aberration causing partial deletion of exon 6 in the *PAH* gene (which in conjunction with a point mutation is causative of phenylketonuria). Data kindly provided by Madhuri Hegde, Ph.D., FACMG, Emory University.

## Probe design

To detect these micro aberrations, highly optimised probes are required. The probes have been selected from OGT's proprietary Oligome™ database — a database of over 25 million in silico optimised oligonucleotide probes. CytoSure arrays utilise 60-mer oligonucleotide probes, which have been shown to offer higher signal-to-noise ratios through increased specificity and sensitivity<sup>2</sup>. To further improve performance, during the design process multiple probes are designed for each of the target regions. These probes are tested in competitive hybridisation experiments and ranked based on technical performance — with only the best performing probes being chosen for the final array design. In addition, further experimental validation using clinical research samples has been performed by Emory University who collaborated in the design and optimisation of CytoSure Disease-Focused Research arrays.

### Cost-effective, customisable analysis

Genetic analysis is typically done in batches to ensure cost-effective processing. However, collecting sufficient samples for rare disorders can significantly delay time to reporting. Through combining informative probes for multiple molecular disorders on a single array, CytoSure Disease-Focused Research arrays allow immediate and accurate processing of diverse samples, thereby reducing both time-to-results and cost.

Each CytoSure Disease-Focused Research array can be modified to create bespoke custom arrays, to suit your specific requirements. The customisation process is made easy with a dedicated project manager from our team of experienced computational biologists assigned to each new project. Additional content, specific to your research interest, can easily be selected from any existing array, our Oligome™ database of pre-optimised probes, or designed de novo to your specifications using our proprietary probe design pipeline. Delivery of your custom array only occurs once you are completely satisfied with the design. A number of array formats are available to suit your design and throughput requirements, with the most typical being 8x60k or 4x180k.

### Easy data interpretation using industry-leading CytoSure Interpret Software

All CytoSure arrays are provided with CytoSure Interpret Software which is a powerful, easy-to-use package for the analysis of aCGH data. Innovative features such as the Accelerate Workflow enable standardised and automated data analysis, including automatic aberration detection and classification. In addition, extensive annotation tracks covering syndromes, genes, exons, CNV and recombination hotspots — each of which link to publically available databases — can be used to provide results in context (Figure 3).




Figure 3: Fully customisable tracks in CytoSure Interpret Software simplify interpretation of aberrations.

### Optimised labelling kits and integrated sample tracking probes

All CytoSure arrays have been validated using CytoSure Genomic DNA Labelling Kits; these labelling kits have been developed and optimised to enable rapid delivery of high-quality results with excellent signal-to-noise ratios. The kits offer faster and simpler DNA labelling and clean-up than alternative enzymatic labelling procedures with improved data quality. Two kit formats are available; the CytoSure Genomic DNA Labelling Kit for 24 samples processed in tubes or the CytoSure HT Genomic DNA Labelling Kit for plate-based processing of up to 96 samples. For best results, CytoSure Genomic DNA Labelling Kits should be used with all CytoSure arrays to give the best derivative log ratio spread (DLRS) values and signal-to-noise ratios ensuring accurate detection of even the smallest aberration.

Increasing numbers of aCGH samples combined with higher-throughput array formats means that it is imperative to track samples throughout the labelling, hybridisation and analysis process to maintain sample identity. CytoSure Sample Tracking Spike-ins are uniquely designed to enable reliable sample tracking and easy identification of sample mix-up using OGT's CytoSure Disease-Focused Research arrays and class-leading CytoSure Interpret Software.

Product	Cat. No.	Format	Genes targeted	Example of diseases covered
CytoSure Autism Research array	700121	4x180k	227	Autism, hearing loss, X-linked intellectual disability.
CytoSure Epilepsy Research array	700112	4x180k	212	Epilepsy, brain malformations, severe combined immunodeficiency.
CytoSure NMD Research array	700117	4x180k	205	Duchenne muscular dystrophy, Becker muscular dystrophy, limb-girdle muscular dystrophy, congenital muscular dystrophy.
CytoSure Cardiomyopathy Research array	700110	4x180k	223	Long QT Syndrome, dilated cardiomyopathy, left ventricular non-compaction, hereditary neuropathies, connective tissue disorders.
CytoSure Eye Disease Research array	700113	4x180k	221	Retinitis pigmentosa, Stargardt disease, macular dystrophy.
CytoSure Hereditary Cancer Research array	700115	4x180k	228	Hereditary cancers, Kabuki syndrome, Lynch syndrome.
CytoSure Metabolic Disorder Research array	700116	4x180k	203	Glycogen storage disease, lysosomal storage disease, mitochondrial disorders (nuclear genes only).
CytoSure Skeletal Dysplasia Research array	700118	4x180k	234	Skeletal dysplasia, limb malformation, severe combined immunodeficiency, osteogenesis imperfecta.
CytoSure Ciliopathy Research array	700111	4x180k	207	Ciliopathies, Stickler syndrome, nephronophthisis, tuberous sclerosis.
CytoSure DMD Research array	Enquire	8x60k	50	Duchenne muscular dystrophy, Becker muscular dystrophy.

### Ordering information

UK +44 (0) 1865 856800

US +1 914 467 5285

contact@ogt.com

ogt.com

Product	Contents	Cat. No.
Disease-Focused Research arrays	See Table 1 for information about each array	Various
CytoSure HT Genomic DNA Labelling Kit	96 reactions, 2 purification plates, dyes, nucleotide mix, random primers, enzyme	500040
CytoSure Sample Tracking Spike-ins A – H	Sample Tracking Probes sufficient for 12 reactions supplied in three aliquots	500066

### References

1. Detection limit of intragenic deletions with targeted array comparative genomic hybridization. BMC Genetics, 2013, 14:116
2. The pitfalls of platform comparison: DNA copy number array technologies assessed. BMC Genomics 2009, 10:588



A Sysmex Group Company

**What binds us,  
makes us.**

Oxford Gene Technology Ltd., Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK

**CytoSure:** This product is provided under an agreement between Agilent Technologies, Inc. and OGT. The manufacture, use, sale or import of this product may be subject to one or more of U.S. patents, pending applications, and corresponding international equivalents, owned by Agilent Technologies, Inc. The purchaser has the non-transferable right to use and consume the product for RESEARCH USE ONLY AND NOT for DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used, for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose. This document and its contents are © Oxford Gene Technology IP Limited – 2021. All rights reserved. OGT™, Oligome™ and CytoSure® are trademarks of Oxford Gene Technology IP Limited.