

SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Hazardous Substances (Safety Data Sheets) Notice 2017. This notice is issued by the Environmental Protection Authority
under sections 75 and 76(1)(b), (f), (g) and (h) of the Hazardous Substances and New Organisms Act 1996

Issuing Date: 10-Aug-23

Revision Date: 10-Aug-23

Revision Number: 3

Section 1: Identification

Product identifier

Product Name DAPI Antifade Counterstain

Product Code(s) CE-DES 500L / CE-DES 1000L
DES 150L / DES 300L / DES 500L / DES 600L / DES 1000L / DES 1200L / DSS 500L /
DFS 500L
RU-DES 150L / RU-DES 300L / RU-DES 500L / RU-DES 600L / RU-DES 1000L / RU-DES
1200L / RU-DSS 500L / RU-DFS 500L

Other means of identification

Recommended use of the chemical and restrictions on use

Recommended use Laboratory chemicals
For professional use only

Uses advised against None known

Details of the supplier of the safety data sheet

Importer

Sysmex New Zealand Limited
Level 3, 103 Carlton Gore Rd
New Market
Auckland 1023, New Zealand
+64-9-630-3554/ 0800797639

Manufacturer

CytoCell Ltd., Oxford Gene Technology
418 Cambridge Science Park, Milton Road,
Cambridge
CB4 0PZ, United Kingdom
T: +44 (0)1223 294048
F: +44 (0)1223 294986
probes@cytoCell.com

E-mail address regulatory@sysmex.co.nz

Emergency telephone number

Emergency telephone For Sysmex Supply Chain support or Product Related Enquiries: +64 9 6303554 /
0800797639 (Mon to Fri – 8.30 am to 5.00 pm)
For any spillage or clean up issues: CHEMCALL 0800 243 622 (24 hours – 365 days)
National Poison Centre 0800 764 766 (0800 POISON)

Section 2: Hazard identification

GHS Classification

Not classified

Label elements

Hazard statements
Not classified

Precautionary Statements - Storage
Store in a well-ventilated place. Keep cool

Precautionary Statements - Disposal
Dispose of contents/container in accordance with local, regional, national, and international regulations as applicable

Other hazards which do not result in classification
No information available.

Section 3: Composition/information on ingredients

Chemical name	CAS No	Weight-%
Glycerol	56-81-5	>=80
Non-hazardous ingredients	Proprietary	Balance

Section 4: First-aid measuresDescription of first aid measures

Inhalation	Remove to fresh air.
Eye contact	Rinse thoroughly with plenty of water, also under the eyelids.
Skin contact	Wash skin with soap and water.
Ingestion	Rinse mouth.

Most important symptoms and effects, both acute and delayed

Symptoms No information available.

Indication of any immediate medical attention and special treatment needed

Note to doctors Treat symptomatically.

Section 5: Fire-fighting measuresSuitable Extinguishing Media

Suitable Extinguishing Media	Dry chemical, CO2, alcohol-resistant foam or water spray.
Unsuitable extinguishing media	None known.

Specific hazards arising from the chemical

Specific hazards arising from the chemical	No information available.
Hazardous combustion products	Carbon oxides.

Special protective actions for fire-fighters

Special protective equipment and precautions for fire-fighters Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

Section 6: Accidental release measuresPersonal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation.
 For emergency responders Use personal protection recommended in Section 8.

Environmental precautions

Environmental precautions See Section 12 for additional Ecological Information. Do not allow to enter into surface water or drains.

Methods and material for containment and cleaning up

Methods for containment Absorb with earth, sand or other non-combustible material and transfer to containers for later disposal. Prevent further leakage or spillage if safe to do so.
 Methods for cleaning up Pick up and transfer to properly labelled containers.

Precautions to prevent secondary hazards

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

Section 7: Handling and storagePrecautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or smoke when using this product. Wear personal protective equipment. Wash hands thoroughly after handling.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep in a dry, cool and well-ventilated place. Protect from moisture.
 Incompatible materials Strong oxidising agents.

Section 8: Exposure controls/personal protectionControl parameters

Exposure Limits

Chemical name	New Zealand	ACGIH TLV	United Kingdom	Australia
Glycerol 56-81-5	TWA: 10 mg/m ³	-	TWA: 10 mg/m ³ STEL: 30 mg/m ³	10 mg/m ³ TWA

Biological occupational exposure limits This product, as supplied, does not contain any hazardous materials with biological limits

established by the region specific regulatory bodies.

Appropriate engineering controls

Engineering controls	Showers Eyewash stations Ventilation systems.
----------------------	---

Individual protection measures, such as personal protective equipment

Eye/face protection	Handling of larger amounts: Wear protective eye glasses for protection against liquid splashes.
Hand protection	To protect the wearer, gloves must be the correct fit and be used properly. Ensure that the breakthrough time of the glove material is not exceeded. Refer to glove supplier for information on breakthrough time for specific gloves.
Skin and body protection	Wear suitable protective clothing.
Respiratory protection	No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.
Environmental exposure controls	No information available.

Section 9: Physical and chemical properties

Information on basic physical and chemical properties

Appearance	Colored liquid
Physical state	Liquid
Colour	Purple
Odour	Odourless
Odour threshold	No information available

Values

		<u>Remarks • Method</u>
pH		No data available
Melting point / freezing point		No data available
Initial boiling point and boiling range	290 °C	
Flash point	176 °C	
Evaporation rate		No data available
Flammability		Not applicable
Flammability Limit in Air		
Upper flammability or explosive limits		Not applicable
Lower flammability or explosive limits		Not applicable
Vapour pressure		No data available
Vapour density		No data available
Relative density		No data available
Water solubility		No data available
Solubility(ies)		No data available
Partition coefficient		No data available
Autoignition temperature		No data available
Decomposition temperature		No data available
Kinematic viscosity		No data available
Dynamic viscosity		No data available
Explosive properties	No information available.	
Oxidising properties	No information available.	

Other information

Softening point	No information available
Molecular weight	No information available
VOC Content (%)	No information available
Liquid Density	No information available
Bulk density	No information available
Particle characteristics	

Section 10: Stability and reactivityReactivity

Reactivity	None under normal use conditions.
------------	-----------------------------------

Chemical stability

Stability	Stable under normal conditions.
-----------	---------------------------------

Explosion data

Sensitivity to mechanical impact	None.
Sensitivity to static discharge	None.

Possibility of hazardous reactions

Possibility of hazardous reactions	None under normal processing.
------------------------------------	-------------------------------

Conditions to avoid

Conditions to avoid	Heat. Protect from direct sunlight. Moisture.
---------------------	---

Incompatible materials

Incompatible materials	Strong oxidising agents.
------------------------	--------------------------

Hazardous decomposition products

Hazardous decomposition products	Carbon oxides.
----------------------------------	----------------

Section 11: Toxicological informationAcute toxicityInformation on likely routes of exposure

Inhalation	Specific test data for the substance or mixture is not available.
Eye contact	Specific test data for the substance or mixture is not available.
Skin contact	Specific test data for the substance or mixture is not available.
Ingestion	Specific test data for the substance or mixture is not available.
Symptoms	No information available.

Acute toxicity

Numerical measures of toxicity
 No information available.

The following values are calculated based on chapter 3.1 of the GHS document:
 ATEmix (oral) 15,750.00 mg/kg

Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Glycerol	= 12600 mg/kg (Rat)	> 10 g/kg (Rabbit)	> 2.75 mg/L (Rat) 4 h

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation No information available.
 Serious eye damage/eye irritation No information available.
 Respiratory or skin sensitisation No information available.
 Germ cell mutagenicity No information available.
 Carcinogenicity No information available.
 Reproductive toxicity No information available.
 STOT - single exposure No information available.
 STOT - repeated exposure No information available.
 Aspiration hazard No information available.

Data used to identify the health effects Refer to Section 16 for Key literature references and sources for data used to compile the SDS.

Section 12: Ecological information

Ecotoxicity

Ecotoxicity The environmental impact of this product has not been fully investigated.

Aquatic ecotoxicity

Chemical name	Algae/aquatic plants	Fish	Crustacea
Glycerol	-	LC50: 51 - 57mL/L (96h, Oncorhynchus mykiss)	-

Terrestrial ecotoxicity There is no data for this product.

Persistence and degradability No information available.

Bioaccumulative potential

Bioaccumulation

Component Information

Chemical name	Partition coefficient
Glycerol	-1.76

Mobility in soil

Mobility in soil No information available.

Other adverse effects

No information available.

Section 13: Disposal considerations

Waste treatment methods

Waste from residues/unused products Not applicable. Not Hazardous.

Contaminated packaging Not applicable. Not Hazardous.

Section 14: Transport information

IATA Not regulated

IMDG Not regulated

Section 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

EPA New Zealand HSNO approval code or group standard Not applicable

National regulations There are no applicable tolerable exposure limits or environmental exposure limits according to the EPA Controls for Hazardous Substances

Certified handlers, tracking and controlled substance license requirements Certified handlers are required for some substances. This includes substances requiring a controlled substance license, and most explosives, vertebrates toxic agents, and certain fumigants. Acutely toxic substances which are a Category 1 or 2, such as pesticides also require Certified handlers. Please check the Health and Safety at Work Act 2015 for further information
Tracking is required for some highly hazardous substances. These substances need to be under the control of an appropriately trained person or appropriately secured. Please check the Health and Safety at Work Act 2015 for further information
Controlled substance licenses are required to possess certain explosives, vertebrate toxic agents and fumigants. See Part 7 of the Health and Safety at Work Regulation 2017 for more information

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

International Inventories

NZIoC	Contact supplier for inventory compliance status.
TSCA	Contact supplier for inventory compliance status.
DSL/NDSL	Contact supplier for inventory compliance status.
EINECS/ELINCS	Contact supplier for inventory compliance status.
ENCS	Contact supplier for inventory compliance status.
IECSC	Contact supplier for inventory compliance status.
KECL	Contact supplier for inventory compliance status.
PICCS	Contact supplier for inventory compliance status.
AICS	Contact supplier for inventory compliance status.

Legend:

- NZIoC - New Zealand Inventory of Chemicals
- TSCA - United States Toxic Substances Control Act Section 8(b) Inventory
- DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List
- EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
- ENCS - Japan Existing and New Chemical Substances
- IECSC - China Inventory of Existing Chemical Substances
- KECL - Korean Existing and Evaluated Chemical Substances
- PICCS - Philippines Inventory of Chemicals and Chemical Substances
- AICS - Australian Inventory of Chemical Substances

Section 16: Other information

Issuing Date	14-Oct-22
Revision Date	14-Oct-22
Revision Note	Initial Release.

Key or legend to abbreviations and acronyms used in the safety data sheet

Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

TWA	TWA (time-weighted average)	STEL	STEL (Short Term Exposure Limit)
Ceiling	Maximum limit value	*	Skin designation
C	Carcinogen		

Key literature references and sources for data used to compile the SDS

- Agency for Toxic Substances and Disease Registry (ATSDR)
- U.S. Environmental Protection Agency ChemView Database
- European Food Safety Authority (EFSA)
- EPA (Environmental Protection Agency)
- Acute Exposure Guideline Level(s) (AEGl(s))
- U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
- U.S. Environmental Protection Agency High Production Volume Chemicals
- Food Research Journal
- Hazardous Substance Database
- International Uniform Chemical Information Database (IUCLID)
- Japan GHS Classification
- Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)
National Library of Medicine's ChemID Plus (NLM CIP)
National Library of Medicine's PubMed database (NLM PUBMED)
National Toxicology Program (NTP)
New Zealand's Chemical Classification and Information Database (CCID)
Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications
Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme
Organisation for Economic Co-operation and Development Screening Information Data Set
World Health Organization

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

End of Safety Data Sheet