



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:  
Regulation (EC) No. 1907/2006 as amended by Commission Regulation (EU) 2020/878 and  
Regulation (EC) No. 1272/2008

Issuing Date 17-Aug-2023    Revision Date 17-Aug-2023    Revision Number 1

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

**Product Code(s)** PCN003  
**Product Name** Mounting Medium  
**Synonyms** None  
**Pure substance/mixture** Mixture

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Recommended use** Laboratory chemicals  
For research use only  
**Uses advised against** None known

### 1.3. Details of the supplier of the safety data sheet

#### Supplier

Cytocell Ltd., Oxford Gene Technology  
418 Cambridge Science Park, Milton Road,  
Cambridge  
CB4 0PZ, United Kingdom  
T: +44 (0)1223 294048  
F: +44 (0)1223 294986  
probes@cytocell.com  
<http://www.ogt.com>

#### For further information, please contact

**E-mail address** probes@cytocell.com

### 1.4. Emergency telephone number

**Emergency telephone** +44 (0) 1223 294048 (Monday - Friday, 9am - 5pm)

**Emergency telephone - §45 - (EC)1272/2008**

**Europe** 112

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]  
This mixture is classified as not hazardous according to regulation (EC) 1272/2008 [CLP].

### 2.2. Label elements

#### Hazard statements

Not classified.

**2.3. Other hazards**

The product does not contain any substance(s) classified as PBT or vPvB.

**Endocrine Disruptor Information** This product does not contain any known or suspected endocrine disruptors.

**SECTION 3: Composition/information on ingredients****3.1 Substances**

Chemical name	Weight-%	REACH registration number	EC No (EU Index No)	Classification according to Regulation (EC) No. 1272/2008 [CLP]	Specific concentration limit (SCL)	M-Factor	M-Factor (long-term)
Glycerol 56-81-5	>80	No data available	200-289-5	[C]	-	-	-

*Classification according to Regulation (EC) No. 1272/2008 [CLP] - Notes*

[C] - Components with occupational exposure limits and/or biological occupational exposure limits requiring monitoring

**Full text of H- and EUH-phrases: see section 16**

**Acute Toxicity Estimate**

If LD50/LC50 data is not available or does not correspond to the classification category, then the appropriate conversion value from CLP Annex I, Table 3.1.2, is used to calculate the acute toxicity estimate (ATEmix) for classifying a mixture based on its components

Chemical name	Oral LD50 mg/kg	Dermal LD50 mg/kg	Inhalation LC50 - 4 hour - dust/mist - mg/L	Inhalation LC50 - 4 hour - vapour - mg/L	Inhalation LC50 - 4 hour - gas - ppm
Glycerol 56-81-5	12600	10010	2.7527	No data available	No data available

This product does not contain candidate substances of very high concern at a concentration  $\geq 0.1\%$  (Regulation (EC) No. 1907/2006 (REACH), Article 59)

**SECTION 4: First aid measures****4.1. Description of first aid measures**

<b>General advice</b>	Get medical attention if irritation or other symptoms occur. Show this safety data sheet to the doctor in attendance.
<b>Inhalation</b>	Remove person to fresh air and keep comfortable for breathing.
<b>Eye contact</b>	Rinse thoroughly with plenty of water, also under the eyelids.
<b>Skin contact</b>	Wash skin with soap and water.
<b>Ingestion</b>	Rinse mouth. Do not induce vomiting without medical advice.

**4.2. Most important symptoms and effects, both acute and delayed**

**Symptoms** No information available.

**Effects of Exposure** No information available.

#### **4.3. Indication of any immediate medical attention and special treatment needed**

**Note to doctors** Treat symptomatically.

### **SECTION 5: Firefighting measures**

#### **5.1. Extinguishing media**

**Suitable Extinguishing Media** Dry chemical, CO2, alcohol-resistant foam or water spray.

**Unsuitable extinguishing media** Do not scatter spilled material with high pressure water streams.

#### **5.2. Special hazards arising from the substance or mixture**

**Specific hazards arising from the chemical** None known based on information supplied.

**Hazardous combustion products** Carbon oxides.

#### **5.3. Advice for firefighters**

**Special protective equipment and precautions for fire-fighters** Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

### **SECTION 6: Accidental release measures**

#### **6.1. Personal precautions, protective equipment and emergency procedures**

**Personal precautions** Ensure adequate ventilation. Use personal protective equipment as required. See section 8 for more information. Avoid contact with skin, eyes or clothing. Avoid breathing vapours or mists. Do not touch or walk through spilled material.

**For emergency responders** Use personal protection recommended in Section 8.

#### **6.2. Environmental precautions**

**Environmental precautions** Do not allow to enter into surface water or drains.

#### **6.3. Methods and material for containment and cleaning up**

**Methods for containment** Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Cover liquid spill with sand, earth or other noncombustible absorbent material. Pick up and transfer to properly labelled containers.

**Prevention of secondary hazards** Clean contaminated objects and areas thoroughly observing environmental regulations.

#### **6.4. Reference to other sections**

**Reference to other sections** See section 8 for more information See section 13 for more information

### **SECTION 7: Handling and storage**

#### **7.1. Precautions for safe handling**

**Advice on safe handling** Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Avoid breathing vapours or mists. Do not eat, drink or smoke when using this product. Wash thoroughly after handling. Wash hands thoroughly after handling. Wear personal protective clothing (see section 8).

**General hygiene considerations** Handle in accordance with good industrial hygiene and safety practice.

### 7.2. Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

**Storage class (TRGS 510)** LGK 10.

### 7.3. Specific end use(s)

**Specific use(s)** The identified uses for this product are detailed in Section 1.2.

## **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

#### Exposure Limits

Chemical name	European Union	Austria	Belgium	Bulgaria	Croatia
Glycerol 56-81-5	-	-	TWA: 10 mg/m <sup>3</sup>	-	TWA: 10 mg/m <sup>3</sup>
Chemical name	Cyprus	Czech Republic	Denmark	Estonia	Finland
Glycerol 56-81-5	-	TWA: 10 mg/m <sup>3</sup> Ceiling: 15 mg/m <sup>3</sup>	-	TWA: 10 mg/m <sup>3</sup>	TWA: 20 mg/m <sup>3</sup>
Chemical name	France	Germany TRGS	Germany DFG	Greece	Hungary
Glycerol 56-81-5	TWA: 10 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup> Peak: 400 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup>	-
Chemical name	Luxembourg	Malta	Netherlands	Norway	Poland
Glycerol 56-81-5	-	-	-	-	TWA: 10 mg/m <sup>3</sup>
Chemical name	Portugal	Romania	Slovakia	Slovenia	Spain
Glycerol 56-81-5	TWA: 10 mg/m <sup>3</sup>	-	TWA: 11 mg/m <sup>3</sup>	TWA: 200 mg/m <sup>3</sup> STEL: 400 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup>
Chemical name	Sweden		Switzerland	United Kingdom	
Glycerol 56-81-5	-		TWA: 50 mg/m <sup>3</sup> STEL: 100 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup> STEL: 30 mg/m <sup>3</sup>	

#### Biological occupational exposure limits

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

#### Derived No Effect Level (DNEL) - Workers

Chemical name	Oral	Dermal	Inhalation
Glycerol 56-81-5	-	-	56 mg/m <sup>3</sup> [5] [6]

#### Notes

[5] Local health effects.  
[6] Long term.

#### Derived No Effect Level (DNEL) - General Public

Chemical name	Oral	Dermal	Inhalation
Glycerol 56-81-5	229 mg/kg bw/day [4] [6]	-	33 mg/m <sup>3</sup> [5] [6]

**Notes**

[4]	Systemic health effects.
[5]	Local health effects.
[6]	Long term.

**Predicted No Effect Concentration (PNEC)**

Chemical name	Freshwater	Freshwater (intermittent release)	Marine water	Marine water (intermittent release)	Air
Glycerol 56-81-5	0.885 mg/L	8.85 mg/L	0.0885 mg/L	-	-

Chemical name	Freshwater sediment	Marine sediment	Sewage treatment	Soil	Food chain
Glycerol 56-81-5	3.3 mg/kg sediment dw	0.33 mg/kg sediment dw	1000 mg/L	0.141 mg/kg soil dw	-

**8.2. Exposure controls****Engineering controls**

Showers  
Eyewash stations  
Ventilation systems.

**Personal protective equipment****Eye/face protection**

Handling of larger amounts: Wear protective eye glasses for protection against liquid splashes. Eye protection must conform to standard EN 166.

**Hand protection**

Wear suitable gloves. Ensure that the breakthrough time of the glove material is not exceeded. Refer to glove supplier for information on breakthrough time for specific gloves. Gloves must conform to standard EN 374.

**Skin and body protection**

No special protective equipment required.

**Respiratory protection**

No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

**General hygiene considerations**

Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls**

No information available.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties****Appearance****Physical state**

Liquid

**Colour**

Colourless

**Odour** Odourless  
**Odour threshold** No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
<b>Melting point / freezing point</b>		No data available
<b>Initial boiling point and boiling range</b>		No data available
<b>Flammability</b>		Not applicable
<b>Flammability Limit in Air</b>		
<b>Upper flammability or explosive limits</b>		Not applicable
<b>Lower flammability or explosive limits</b>		Not applicable
<b>Flash point</b>		No data available
<b>Autoignition temperature</b>		No data available
<b>Decomposition temperature</b>		No data available
<b>pH</b>		No data available
<b>pH (as aqueous solution)</b>		No data available
<b>Kinematic viscosity</b>		No data available
<b>Dynamic viscosity</b>		No data available
<b>Water solubility</b>	Soluble in water	
<b>Solubility(ies)</b>		No data available
<b>Partition coefficient</b>		No data available
<b>Vapour pressure</b>		No data available
<b>Relative density</b>		No data available
<b>Bulk density</b>		No data available
<b>Liquid Density</b>		No data available
<b>Relative vapour density</b>		No data available
<b>Particle characteristics</b>		
<b>Particle Size</b>		No data available
<b>Particle Size Distribution</b>		No data available

## 9.2. Other information

9.2.1. Information with regards to physical hazard classes  
 Not applicable

9.2.2. Other safety characteristics  
 No information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

**Reactivity** None under normal use conditions.

### 10.2. Chemical stability

**Stability** Stable under normal conditions.

### Explosion data

**Sensitivity to mechanical impact** None.

**Sensitivity to static discharge** None.

### 10.3. Possibility of hazardous reactions

**Possibility of hazardous reactions** None under normal processing.

### 10.4. Conditions to avoid

**Conditions to avoid** Heat. Moisture. Protect from direct sunlight.

**10.5. Incompatible materials**

**Incompatible materials** Strong oxidising agents.

**10.6. Hazardous decomposition products**

**Hazardous decomposition products** None known based on information supplied.

**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008****Information on likely routes of exposure****Product Information**

<b>Inhalation</b>	Specific test data for the substance or mixture is not available.
<b>Eye contact</b>	Specific test data for the substance or mixture is not available.
<b>Skin contact</b>	Non-irritating to the skin.
<b>Ingestion</b>	Specific test data for the substance or mixture is not available.

**Symptoms related to the physical, chemical and toxicological characteristics**

**Symptoms** No information available.

**Acute toxicity****Numerical measures of toxicity**

Based on available data, the classification criteria are not met.

**Component Information**

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Glycerol	= 12600 mg/kg ( Rat )	> 10 g/kg ( Rabbit )	> 2.75 mg/L ( Rat ) 4 h

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

<b>Skin corrosion/irritation</b>	Based on available data, the classification criteria are not met.
<b>Serious eye damage/eye irritation</b>	Based on available data, the classification criteria are not met.
<b>Respiratory or skin sensitisation</b>	Based on available data, the classification criteria are not met.
<b>Germ cell mutagenicity</b>	Based on available data, the classification criteria are not met.
<b>Carcinogenicity</b>	Based on available data, the classification criteria are not met.
<b>Reproductive toxicity</b>	Based on available data, the classification criteria are not met.
<b>STOT - single exposure</b>	Based on available data, the classification criteria are not met.
<b>STOT - repeated exposure</b>	Based on available data, the classification criteria are not met.
<b>Aspiration hazard</b>	Based on available data, the classification criteria are not met.

**11.2. Information on other hazards****11.2.1. Endocrine disrupting properties**

**Endocrine disrupting properties** This product does not contain any known or suspected endocrine disruptors.

**11.2.2. Other information**

**Other adverse effects** No information available.

**SECTION 12: Ecological information****12.1. Toxicity**

**Ecotoxicity** Not considered to be harmful to aquatic life.

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Glycerol 56-81-5	-	LC50: 51 - 57mL/L (96h, Oncorhynchus mykiss)	-	-

**12.2. Persistence and degradability**

**Persistence and degradability** The product is substantially biodegradable.

**12.3. Bioaccumulative potential****Bioaccumulation****Component Information**

Chemical name	Partition coefficient
Glycerol	-1.75

**12.4. Mobility in soil**

**Mobility in soil** Soluble in water. The product is predicted to have high mobility in soil.

**Mobility** Soluble in water.

**12.5. Results of PBT and vPvB assessment**

**PBT and vPvB assessment** The product does not contain any substance(s) classified as PBT or vPvB.

Chemical name	PBT and vPvB assessment
Glycerol 56-81-5	The substance is not PBT / vPvB

**12.6. Endocrine disrupting properties**

**Endocrine disrupting properties** This product does not contain any known or suspected endocrine disruptors.

**12.7. Other adverse effects**

**Other adverse effects** No information available.

**SECTION 13: Disposal considerations**



**13.1. Waste treatment methods**

<b>Waste from residues/unused products</b>	Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.
<b>Contaminated packaging</b>	Do not reuse empty containers.
<b>Waste codes / waste designations according to EWC / AVV</b>	According to the European Waste Catalogue, Waste Codes are not product specific, but application specific. Waste codes should be assigned by the user based on the application for which the product was used.

**SECTION 14: Transport information**

<b>IMDG</b>	Not regulated
<b>14.1 UN number or ID number</b>	Not regulated
<b>14.2 UN proper shipping name</b>	Not regulated
<b>14.3 Transport hazard class(es)</b>	Not regulated
<b>14.4 Packing group</b>	Not applicable
<b>14.5 Environmental hazards</b>	Not applicable
<b>14.6 Special Precautions for Users</b>	
<b>Special Provisions</b>	None
<b>14.7 Maritime transport in bulk according to IMO instruments</b>	No information available

<b>RID</b>	Not regulated
<b>14.1 UN number</b>	Not regulated
<b>14.2 UN proper shipping name</b>	Not regulated
<b>14.3 Transport hazard class(es)</b>	Not regulated
<b>14.4 Packing group</b>	Not applicable
<b>14.5 Environmental hazards</b>	Not applicable
<b>14.6 Special Precautions for Users</b>	
<b>Special Provisions</b>	None

<b>ADR</b>	Not regulated
<b>14.1 UN number or ID number</b>	Not regulated
<b>14.2 UN proper shipping name</b>	Not regulated
<b>14.3 Transport hazard class(es)</b>	Not regulated
<b>14.4 Packing group</b>	Not applicable
<b>14.5 Environmental hazards</b>	Not applicable
<b>14.6 Special Precautions for Users</b>	
<b>Special Provisions</b>	None

<b>IATA</b>	Not regulated
<b>14.1 UN number or ID number</b>	Not regulated
<b>14.2 UN proper shipping name</b>	Not regulated
<b>14.3 Transport hazard class(es)</b>	Not regulated
<b>14.4 Packing group</b>	Not applicable
<b>14.5 Environmental hazards</b>	Not applicable
<b>14.6 Special Precautions for Users</b>	
<b>Special Provisions</b>	None
<b>Note:</b>	None

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

**European Union**

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work.

**Authorisations and/or restrictions on use:**

This product does not contain substances subject to authorisation (Regulation (EC) No. 1907/2006 (REACH), Annex XIV) This product does not contain substances subject to restriction (Regulation (EC) No. 1907/2006 (REACH), Annex XVII)

**Persistent Organic Pollutants**

Not applicable

**Ozone-depleting substances (ODS) regulation (EC) 1005/2009**

Not applicable

**International Inventories**

Contact supplier for inventory compliance status

**15.2. Chemical safety assessment****Chemical Safety Report**

No information available

**SECTION 16: Other information****Key or legend to abbreviations and acronyms used in the safety data sheet****Legend**

ATE: Acute Toxicity Estimate

SVHC: Substances of Very High Concern for Authorisation:

PBT: Persistent, Bioaccumulative, and Toxic (PBT) Substances

vPvB: Very Persistent and very Bioaccumulative (vPvB) Substances

**Legend Section 8: Exposure controls/personal protection**

TWA TWA (time-weighted average)

STEL

STEL (Short Term Exposure Limit)

Ceiling Maximum limit value

\*

Skin designation

SCBA Self-contained breathing apparatus

Classification procedure	
Classification according to Regulation (EC) No. 1272/2008 [CLP]	Method Used
Acute oral toxicity	Calculation method
Acute dermal toxicity	Calculation method
Acute inhalation toxicity - gas	Calculation method
Acute inhalation toxicity - vapour	Calculation method
Acute inhalation toxicity - dust/mist	Calculation method
Skin corrosion/irritation	Calculation method
Serious eye damage/eye irritation	Calculation method
Respiratory sensitisation	Calculation method
Skin sensitisation	Calculation method
Mutagenicity	Calculation method
Carcinogenicity	Calculation method

Reproductive toxicity	Calculation method
STOT - single exposure	Calculation method
STOT - repeated exposure	Calculation method
Acute aquatic toxicity	Calculation method
Chronic aquatic toxicity	Calculation method
Aspiration hazard	Calculation method
Ozone	Calculation method

#### Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)  
 U.S. Environmental Protection Agency ChemView Database  
 European Food Safety Authority (EFSA)  
 European Chemicals Agency (ECHA) Committee for Risk Assessment (ECHA\_RAC)  
 European Chemicals Agency (ECHA) (ECHA\_API)  
 EPA (Environmental Protection Agency)  
 Acute Exposure Guideline Level(s) (AEGl(s))  
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act  
 U.S. Environmental Protection Agency High Production Volume Chemicals  
 Food Research Journal  
 Hazardous Substance Database  
 International Uniform Chemical Information Database (IUCLID)  
 Japan GHS Classification  
 Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)  
 NIOSH (National Institute for Occupational Safety and Health)  
 National Library of Medicine's ChemID Plus (NLM CIP)  
 National Library of Medicine's PubMed database (NLM PUBMED)  
 National Toxicology Program (NTP)  
 New Zealand's Chemical Classification and Information Database (CCID)  
 Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications  
 Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme  
 Organisation for Economic Co-operation and Development Screening Information Data Set  
 World Health Organization

**Issuing Date** 17-Aug-2023

**Revision Date** 17-Aug-2023

**Revision Note** Initial Release.

**This safety data sheet complies with the requirements of Commission Regulation (EU) 2020/878 of 18 June 2020 amending Regulation (EC) No. 1907/2006**

#### Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

**End of Safety Data Sheet**