



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Hazardous Substances (Safety Data Sheets) Notice 2017. This notice is issued by the Environmental Protection Authority
under sections 75 and 76(1)(b), (f), (g) and (h) of the Hazardous Substances and New Organisms Act 1996

Issuing Date: 14-Oct-22 Revision Date: 14-Oct-22 Revision Number: 1

Section 1: Identification

Product identifier

Product Name ENZYME SOLUTION

Product Code(s) LPS100B

Other means of identification

Recommended use of the chemical and restrictions on use

Recommended use Laboratory chemicals

For professional use only

Uses advised against None known

Details of the supplier of the safety data sheet

<u>Importer</u> <u>Manufacturer</u>

Sysmex New Zealand Limited Cytocell Ltd., Oxford Gene Technology Level 3, 103 Carlton Gore Rd 418 Cambridge Science Park, Milton Road,

New Market Cambridge

Auckland 1023, New Zealand CB4 0PZ, United Kingdom +64-9-630-3554/0800797639 T: +44 (0)1223 294048

F: +44 (0)1223 294986 probes@cytocell.com

E-mail address regulatory@sysmex.co.nz

Emergency telephone number

Emergency telephone For Sysmex Supply Chain support or Product Related Enquiries: +64 9 6303554 /

0800797639 (Mon to Fri – 8.30 am to 5.00 pm)

For any spillage or clean up issues: CHEMCALL 0800 243 622 (24 hours – 365 days)

National Poison Centre 0800 764 766 (0800 POISON)

Section 2: Hazard identification

GHS Classification

Respiratory sensitisation Category 1

SDSLPS100B NZL V001 Page 1 / 9

Label elements



Signal word

Danger

Hazard statements

May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary Statements - Prevention

Avoid breathing dust/fume/gas/mist/vapours/spray In case of inadequate ventilation wear respiratory protection

Precautionary Statements - Response

IF INHALED: Remove person to fresh air and keep comfortable for breathing If experiencing respiratory symptoms: Call a POISON CENTER or doctor

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Other hazards which do not result in classification

Causes mild skin irritation.

Section 3: Composition/information on ingredients

Chemical name	CAS No	Weight-%
Pepsin A	9001-75-6	0.1 - 5
Non-hazardous ingredients	Proprietary	Balance

Section 4: First-aid measures

Description of first aid measures

General advice Show this safety data sheet to the doctor in attendance

Inhalation May cause allergic respiratory reaction. If breathing has stopped, give artificial respiration.

Get medical attention immediately. Remove to fresh air. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Get immediate medical advice/attention.

Eye contact Rinse thoroughly with plenty of water, also under the eyelids.

Skin contact Wash with soap and water. May cause an allergic skin reaction. In the case of skin irritation

or allergic reactions see a doctor.

Ingestion May produce an allergic reaction. Do NOT induce vomiting. Rinse mouth. Never give

anything by mouth to an unconscious person. Get immediate medical advice/attention.

Self-protection of the first aider Ensure that medical personnel are aware of the material(s) involved, take precautions to

protect themselves and prevent spread of contamination. Avoid contact with skin, eyes or clothing. Avoid direct contact with skin. Use barrier to give mouth-to-mouth resuscitation. Use personal protective equipment as required. See section 8 for more information.

SDSLPS100B NZL V001 Page 2 / 9

Most important symptoms and effects, both acute and delayed

Symptoms May cause allergy or asthma symptoms or breathing difficulties if inhaled. Coughing and/ or

wheezing. Itching. Rashes. Hives. Prolonged contact may cause redness and irritation.

Indication of any immediate medical attention and special treatment needed

Note to doctorsMay cause sensitisation in susceptible persons. Treat symptomatically.

Section 5: Fire-fighting measures

Suitable Extinguishing Media

Suitable Extinguishing Media Dry chemical, CO2, alcohol-resistant foam or water spray.

Unsuitable extinguishing media Do not scatter spilled material with high pressure water streams.

Specific hazards arising from the chemical

Specific hazards arising from the

chemical

Product is or contains a sensitiser. May cause sensitisation by inhalation and skin contact.

Hazardous combustion products Carbon oxides.

Special protective actions for fire-fighters

Special protective equipment and precautions for fire-fighters

Firefighters should wear self-contained breathing apparatus and full firefighting turnout

gear. Use personal protection equipment.

Section 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures

contact with skin, eyes or clothing. Ensure adequate ventilation. Evacuate personnel to safe

areas. Keep people away from and upwind of spill/leak.

Other information Refer to protective measures listed in Sections 7 and 8.

For emergency responders

Use personal protection recommended in Section 8.

Environmental precautions

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not allow to enter into surface water

or drains. Inform authorities in the event of product spillage to water courses or sewage

systems.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so.

Methods for cleaning up Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder,

sawdust). Pick up and transfer to properly labelled containers.

Precautions to prevent secondary hazards

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

SDSLPS100B NZL V001 Page 3 / 9

Section 7: Handling and storage

Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Avoid contact with

skin, eyes or clothing. Ensure adequate ventilation. Provide extract ventilation to points where emissions occur. In case of insufficient ventilation, wear suitable respiratory equipment. Do not eat, drink or smoke when using this product. Remove contaminated clothing and shoes. Wear personal protective equipment. Wash hands thoroughly after

handling.

General hygiene considerations Avoid contact with skin, eyes or clothing. Wear suitable gloves and eye/face protection. Do

not eat, drink or smoke when using this product. Remove and wash contaminated clothing

and gloves, including the inside, before re-use.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place. Store locked up.

Incompatible materials Strong acids. Strong bases. Alcohols.

Section 8: Exposure controls/personal protection

Control parameters

Exposure Limits This product, as supplied, does not contain any hazardous materials with occupational

exposure limits established by the region specific regulatory bodies.

Biological occupational exposure limits This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

Appropriate engineering controls

Engineering controls Showers

Eyewash stations Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection Handling of larger amounts: Wear protective eye glasses for protection against liquid

splashes. Wear safety glasses with side shields (or goggles).

Hand protection To protect the wearer, gloves must be the correct fit and be used properly. Ensure that the

breakthrough time of the glove material is not exceeded. Refer to glove supplier for

information on breakthrough time for specific gloves. Wear suitable gloves.

Skin and body protection Wear suitable protective clothing.

Respiratory protection In case of insufficient ventilation, wear suitable respiratory equipment.

Environmental exposure controls No information available.

SDSLPS100B NZL V001 Page 4 / 9

Section 9: Physical and chemical properties

Information on basic physical and chemical properties

Appearance

Physical state Liquid
Colour Colourless
Odour Odourless

Odour threshold No information available

<u>Values</u> <u>Remarks • Method</u>

pHNo data availableMelting point / freezing pointNo data availableInitial boiling point and boilingNo data available

range

Flash pointNo data availableEvaporation rateNo data availableFlammabilityNot flammable

Flammability Limit in Air

Upper flammability or explosive Not applicable

limits

Lower flammability or explosive Not applicable

limits

Vapour pressureNo data availableVapour densityNo data availableRelative densityNo data available

Water solubility
Solubility(ies)
No data available
Partition coefficient
No data available
Autoignition temperature
No data available
Decomposition temperature
No data available
Kinematic viscosity
No data available

No information available.

Dynamic viscosity
Explosive properties No information available.

Other information

Oxidising properties

Softening point
Molecular weight
VOC Content (%)
Liquid Density
No information available

Particle characteristics

Section 10: Stability and reactivity

Reactivity

Reactivity None under normal use conditions.

Chemical stability

Stability Stable under normal conditions.

Explosion data

Sensitivity to mechanical impact None.
Sensitivity to static discharge None.

SDSLPS100B NZL V001 Page 5 / 9

No data available

Possibility of hazardous reactions

Possibility of hazardous reactions None under normal processing.

Conditions to avoid

Conditions to avoid Protect from direct sunlight. Extremely high or low temperatures.

Incompatible materials

Incompatible materials Strong acids. Strong bases. Alcohols.

Hazardous decomposition products

Hazardous decomposition products Carbon oxides.

Section 11: Toxicological information

Acute toxicity

Information on likely routes of exposure

Product Information

Inhalation May cause sensitisation in susceptible persons. (based on components).

Eye contact May cause slight irritation. (based on components).

Skin contact Repeated or prolonged skin contact may cause allergic reactions with susceptible persons.

(based on components). Causes mild skin irritation.

Ingestion Specific test data for the substance or mixture is not available.

Symptoms Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling

of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain, or flushing. Coughing and/ or wheezing. Prolonged contact may cause redness and irritation.

Acute toxicity

Numerical measures of toxicity

No information available.

The following values are calculated based on chapter 3.1 of the GHS document:

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation Classification based on data available for ingredients. May cause skin irritation.

Serious eye damage/eye irritation No information available.

Respiratory or skin sensitisation May cause sensitisation by inhalation.

Germ cell mutagenicity
Carcinogenicity
Reproductive toxicity
STOT - single exposure
STOT - repeated exposure
Aspiration hazard
No information available.
No information available.
No information available.
No information available.

SDSLPS100B NZL V001 Page 6 / 9

Data used to identify the health effects

Refer to Section 16 for Key literature references and sources for data used to compile the SDS.

Section 12: Ecological information

Ecotoxicity

Ecotoxicity Not considered to be harmful to aquatic life.

Aquatic ecotoxicity

Terrestrial ecotoxictyThere is no data for this product.

Persistence and degradability No information available.

Bioaccumulative potential

Bioaccumulation No information available.

Mobility in soil

Mobility in soil No information available.

Other adverse effects

No information available.

Section 13: Disposal considerations

Waste treatment methods

Waste from residues/unused products

Dispose of product in packaging in a way that is consistent with the EPA Consolidation 30 April 2021 of the Hazardous Substances (Disposal) Notice 2017 and the Act. Treat the substance using a method that changes the characteristics or composition of the substance so that the substance is no longer a hazardous substance; or export the substance from New Zealand as waste. Substances which are hazardous to human health or corrosive to metals – may be discharged into the environment if a tolerable exposure limit has been set for the substance (or a component of that substance); and the discharge does not, after reasonable mixing, result in the concentration of the substance in an environmental medium exceeding the tolerable exposure limit. If there is no tolerable exposure limit for the substance, then it may only be discharged into the environment if the substance is very rapidly converted to substances that are not hazardous substances.

Contaminated packaging Not applicable. Not Hazardous.

Section 14: Transport information

IATA Not regulated

IMDG Not regulated

SDSLPS100B NZL V001 Page 7 / 9

Section 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

EPA New Zealand HSNO approval

code or group standard

HSR002596 - Laboratory Chemicals and Reagent Kits

National regulations

There are no applicable tolerable exposure limits or environmental exposure limits

according to the EPA Controls for Hazardous Substances

Certified handlers, tracking and controlled substance license requirements

Certified handlers are required for some substances. This includes substances requiring a controlled substance license, and most explosives, vertebrates toxic agents, and certain fumigants. Acutely toxic substances which are a Category 1 or 2, such as pesticides also require Certified handlers. Please check the Health and Safety at Work Act 2015 for further information

Tracking is required for some highly hazardous substances. These substances need to be under the control of an appropriately trained person or appropriately secured. Please check

the Health and Safety at Work Act 2015 for further information

Controlled substance licenses are required to possess certain explosives, vertebrate toxic agents and fumigants. See Part 7 of the Health and Safety at Work Regulation 2017 for

more information

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

International Inventories

NZIoC Contact supplier for inventory compliance status. **TSCA** Contact supplier for inventory compliance status. Contact supplier for inventory compliance status. **DSL/NDSL** Contact supplier for inventory compliance status. **EINECS/ELINCS** Contact supplier for inventory compliance status. **ENCS IECSC** Contact supplier for inventory compliance status. Contact supplier for inventory compliance status. KECL **PICCS** Contact supplier for inventory compliance status. **AICS** Contact supplier for inventory compliance status.

Legend:

NZIoC - New Zealand Inventory of Chemicals

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

Section 16: Other information

Issuing Date 14-Oct-22

Revision Date 14-Oct-22

Revision Note Initial Release.

SDSLPS100B NZL V001 Page 8 / 9

Key or legend to abbreviations and acronyms used in the safety data sheet

Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

TWA TWA (time-weighted average) STEL STEL (Short Term Exposure Limit)

Ceiling Maximum limit value * Skin designation

C Carcinogen

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR) U.S. Environmental Protection Agency ChemView Database

European Food Safety Authority (EFSA)

EPA (Environmental Protection Agency)

Acute Exposure Guideline Level(s) (AEGL(s))

U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act

U.S. Environmental Protection Agency High Production Volume Chemicals

Food Research Journal

Hazardous Substance Database

International Uniform Chemical Information Database (IUCLID)

Japan GHS Classification

Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)

National Library of Medicine's ChemID Plus (NLM CIP)

National Library of Medicine's PubMed database (NLM PUBMED)

National Toxicology Program (NTP)

New Zealand's Chemical Classification and Information Database (CCID)

Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications

Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme

Organisation for Economic Co-operation and Development Screening Information Data Set

World Health Organization

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

End of Safety Data Sheet

SDSLPS100B NZL V001 Page 9 / 9