



# SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Hazardous Substances (Safety Data Sheets) Notice 2017. This notice is issued by the Environmental Protection Authority
under sections 75 and 76(1)(b), (f), (g) and (h) of the Hazardous Substances and New Organisms Ac 1996

Issuing Date: 14-Oct-22 Revision Date: 14-Oct-22 Revision Number: 1

### Section 1: Identification

### **Product identifier**

Product Name HEAT PRETREATMENT SOLUTION PH 7

Product Code(s) LPS100A

Other means of identification

Recommended use of the chemical and restrictions on use

Recommended use Laboratory chemicals

For professional use only

Uses advised against None known

#### Details of the supplier of the safety data sheet

<u>Importer</u> <u>Manufacturer</u>

Sysmex New Zealand Limited Cytocell Ltd., Oxford Gene Technology Level 3, 103 Carlton Gore Rd 418 Cambridge Science Park, Milton Road,

New Market Cambridge

F: +44 (0)1223 294986 probes@cytocell.com

E-mail address regulatory@sysmex.co.nz

**Emergency telephone number** 

Emergency telephone For Sysmex Supply Chain support or Product Related Enquiries: +64 9 6303554 /

0800797639 (Mon to Fri – 8.30 am to 5.00 pm)

For any spillage or clean up issues: CHEMCALL 0800 243 622 (24 hours – 365 days)

National Poison Centre 0800 764 766 (0800 POISON)

### Section 2: Hazard identification

### **GHS Classification**

Not classified

SDSLPS100A NZL V001 Page 1 / 8

Revision Date: 14-Oct-22

#### **Label elements**

#### **Hazard statements**

Not classified

### **Precautionary Statements - Storage**

Store in a well-ventilated place. Keep cool

#### **Precautionary Statements - Disposal**

Dispose of contents/container in accordance with local, regional, national, and international regulations as applicable

#### Other hazards which do not result in classification

No information available.

# Section 3: Composition/information on ingredients

**Component Information** The product contains no substances which at their given concentration, are considered to

be hazardous to health

Non-hazardous ingredients Proprietary Balance

### Section 4: First-aid measures

### Description of first aid measures

**Inhalation** Remove to fresh air.

**Eye contact** Rinse thoroughly with plenty of water, also under the eyelids.

**Skin contact** Wash skin with soap and water.

**Ingestion** Rinse mouth.

### Most important symptoms and effects, both acute and delayed

**Symptoms** No information available.

### Indication of any immediate medical attention and special treatment needed

# Section 5: Fire-fighting measures

### Suitable Extinguishing Media

**Suitable Extinguishing Media** Dry chemical, CO2, alcohol-resistant foam or water spray.

Unsuitable extinguishing media None known.

### Specific hazards arising from the chemical

Specific hazards arising from the

chemical

No information available.

**Hazardous combustion products** Carbon oxides.

SDSLPS100A NZL V001 Page 2 / 8

#### Special protective actions for fire-fighters

Special protective equipment and precautions for fire-fighters

Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

Revision Date: 14-Oct-22

### Section 6: Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

**Personal precautions** Ensure adequate ventilation.

**Environmental precautions** 

**Environmental precautions**See Section 12 for additional Ecological Information. Do not allow to enter into surface

water or drains.

#### Methods and material for containment and cleaning up

**Methods for containment**Absorb with earth, sand or other non-combustible material and transfer to containers for

later disposal. Prevent further leakage or spillage if safe to do so.

**Methods for cleaning up** Pick up and transfer to properly labelled containers.

#### Precautions to prevent secondary hazards

**Prevention of secondary hazards** Clean contaminated objects and areas thoroughly observing environmental regulations.

### Section 7: Handling and storage

#### Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or

smoke when using this product. Wear personal protective equipment. Wash hands

thoroughly after handling.

#### Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep in a dry, cool and well-ventilated place.

**Incompatible materials** Strong bases. Strong acids.

# Section 8: Exposure controls/personal protection

#### **Control parameters**

#### **Exposure Limits**

**Biological occupational exposure limits** This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies.

### Appropriate engineering controls

Engineering controls Showers

Eyewash stations Ventilation systems.

SDSLPS100A NZL V001 Page 3 / 8

### Individual protection measures, such as personal protective equipment

Eye/face protection Handling of larger amounts: Wear protective eye glasses for protection against liquid

splashes.

Hand protection To protect the wearer, gloves must be the correct fit and be used properly. Ensure that the

breakthrough time of the glove material is not exceeded. Refer to glove supplier for

Revision Date: 14-Oct-22

information on breakthrough time for specific gloves.

**Skin and body protection** Wear suitable protective clothing.

exceeded or irritation is experienced, ventilation and evacuation may be required.

**Environmental exposure controls** No information available.

# Section 9: Physical and chemical properties

### Information on basic physical and chemical properties

Appearance Colourless liquid

Physical state
Colour
Colourless
Odour
Characteristic

Odour threshold No information available

Values
pH 6 - 8

Melting point / freezing point

No data available
Initial boiling point and boiling

No data available

range
Flash point
Evaporation rate
No data available

Flammability Limit in Air

Upper flammability or explosive Not applicable

limits

Lower flammability or explosive Not applicable

limits

Vapour pressure No data available Vapour density No data available Relative density No data available Water solubility No data available No data available Solubility(ies) Partition coefficient No data available **Autoignition temperature** No data available **Decomposition temperature** No data available Kinematic viscosity No data available **Dvnamic viscosity** No data available

**Explosive properties**No information available. **Oxidising properties**No information available.

**Other information** 

Softening point
Molecular weight
VOC Content (%)
Liquid Density
Bulk density
No information available
No information available
No information available
No information available

Particle characteristics

SDSLPS100A NZL V001 Page 4 / 8

Revision Date: 14-Oct-22

# Section 10: Stability and reactivity

Reactivity

**Reactivity** None under normal use conditions.

**Chemical stability** 

**Stability** Stable under normal conditions.

**Explosion data** 

Sensitivity to mechanical impact None.

Sensitivity to static discharge None.

Possibility of hazardous reactions

Possibility of hazardous reactions 
None under normal processing.

**Conditions to avoid** 

**Conditions to avoid** Protect from direct sunlight. Extremely high or low temperatures.

**Incompatible materials** 

Incompatible materials Strong bases. Strong acids.

**Hazardous decomposition products** 

Hazardous decomposition products Carbon oxides.

# Section 11: Toxicological information

#### **Acute toxicity**

### Information on likely routes of exposure

**Inhalation** Specific test data for the substance or mixture is not available.

**Eye contact** Specific test data for the substance or mixture is not available.

**Skin contact** Specific test data for the substance or mixture is not available.

**Ingestion** Specific test data for the substance or mixture is not available.

**Symptoms** No information available.

**Acute toxicity** 

**Numerical measures of toxicity** 

No information available.

The following values are calculated based on chapter 3.1 of the GHS document:

**ATEmix (oral)** 15,750.00 mg/kg

SDSLPS100A NZL V001 Page 5 / 8

#### **Component Information**

#### Delayed and immediate effects as well as chronic effects from short and long-term exposure

**Skin corrosion/irritation**No information available.

**Serious eye damage/eye irritation** No information available.

**Respiratory or skin sensitisation** No information available.

Germ cell mutagenicity No information available.

**Carcinogenicity** No information available.

**Reproductive toxicity** No information available.

**STOT - single exposure**STOT - repeated exposure
No information available.
No information available.

**Aspiration hazard** No information available.

Data used to identify the health

effects

Refer to Section 16 for Key literature references and sources for data used to compile the

Revision Date: 14-Oct-22

SDS.

# Section 12: Ecological information

#### **Ecotoxicity**

**Ecotoxicity** The environmental impact of this product has not been fully investigated.

**Aquatic ecotoxicity** 

**Terrestrial ecotoxicty** There is no data for this product.

Persistence and degradability No information available.

**Bioaccumulative potential** 

**Bioaccumulation** 

**Component Information** 

**Mobility in soil** 

**Mobility in soil** No information available.

Other adverse effects

No information available.

SDSLPS100A NZL V001 Page 6 / 8

### Section 13: Disposal considerations

#### Waste treatment methods

Waste from residues/unused

products

Not applicable. Not Hazardous.

Contaminated packaging Not applicable. Not Hazardous.

# Section 14: Transport information

**IATA** Not regulated

**IMDG** Not regulated

# Section 15: Regulatory information

#### Safety, health and environmental regulations/legislation specific for the substance or mixture

**EPA New Zealand HSNO approval** 

code or group standard

Not applicable

**National regulations** There are no applicable tolerable exposure limits or environmental exposure limits

according to the EPA Controls for Hazardous Substances

Certified handlers, tracking and controlled substance license requirements

Certified handlers are required for some substances. This includes substances requiring a controlled substance license, and most explosives, vertebrates toxic agents, and certain fumigants. Acutely toxic substances which are a Category 1 or 2, such as pesticides also require Certified handlers. Please check the Health and Safety at Work Act 2015 for further information

Revision Date: 14-Oct-22

Tracking is required for some highly hazardous substances. These substances need to be under the control of an appropriately trained person or appropriately secured. Please check

the Health and Safety at Work Act 2015 for further information

Controlled substance licenses are required to possess certain explosives, vertebrate toxic agents and fumigants. See Part 7 of the Health and Safety at Work Regulation 2017 for more information

#### **International Regulations**

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

**International Inventories** 

**NZIoC** Contact supplier for inventory compliance status. **TSCA** Contact supplier for inventory compliance status. DSL/NDSL Contact supplier for inventory compliance status. **EINECS/ELINCS** Contact supplier for inventory compliance status. Contact supplier for inventory compliance status. **ENCS IECSC** Contact supplier for inventory compliance status. **KECL** Contact supplier for inventory compliance status. **PICCS** Contact supplier for inventory compliance status. **AICS** Contact supplier for inventory compliance status.

SDSLPS100A NZL V001 Page 7 / 8

Revision Date: 14-Oct-22

#### Legend:

**NZIoC** - New Zealand Inventory of Chemicals

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

**IECSC** - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

PICCS - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

### Section 16: Other information

Issuing Date 14-Oct-22

Revision Date 14-Oct-22

Revision Note Initial Release.

#### Key or legend to abbreviations and acronyms used in the safety data sheet

Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

TWA TWA (time-weighted average) STEL STEL (Short Term Exposure Limit)

Ceiling Maximum limit value \* Skin designation

C Carcinogen

#### Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)

U.S. Environmental Protection Agency ChemView Database

European Food Safety Authority (EFSA)

EPA (Environmental Protection Agency)

Acute Exposure Guideline Level(s) (AEGL(s))

U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act

U.S. Environmental Protection Agency High Production Volume Chemicals

Food Research Journal

Hazardous Substance Database

International Uniform Chemical Information Database (IUCLID)

Japan GHS Classification

Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)

NIOSH (National Institute for Occupational Safety and Health)

National Library of Medicine's ChemID Plus (NLM CIP)

National Library of Medicine's PubMed database (NLM PUBMED)

National Toxicology Program (NTP)

New Zealand's Chemical Classification and Information Database (CCID)

Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications

Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme

Organisation for Economic Co-operation and Development Screening Information Data Set

World Health Organization

#### **Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet** 

\_\_\_\_\_

SDSLPS100A NZL V001 Page 8 / 8