

SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
SafeWork Australia Approved Code of Practice about the preparation of safety data sheets for
hazardous chemicals (July 2018), which is an approved code of practice under section 274 of the Work
Health and Safety Act

Issuing Date 11-Jan-23

Revision Date 11-Jan-23

Revision Number 3

Section 1: Identification

Product identifier

Product Name CytoCell and myProbes Liquid FISH Probes

Product Code(s) CE-LP* *** / LP* *** / RU-LP* *** / MP****

Other means of identification

Pure substance/mixture Mixture

Recommended use of the chemical and restrictions on use

Recommended use Laboratory chemicals
For professional use only

Uses advised against None known

Details of manufacturer or importer

Importer

Sysmex Australia Pty Ltd
Suite 3, Level 5
15 Talavera Rd
Macquarie Park
NSW 2113
Telephone no. +61 2 9016 3040

Manufacturer

CytoCell Ltd., Oxford Gene Technology
418 Cambridge Science Park, Milton Road,
Cambridge
CB4 0PZ, United Kingdom
T: +44 (0)1223 294048
F: +44 (0)1223 294986
probes@cytoCell.com

For further information, please contact

Contact Point Product Safety Department

Emergency telephone number

Emergency telephone number For medical advice (English): 13 11 26 (Poisons Information Centre)

Section 2: Hazard(s) identification

Classified as a hazardous substance in accordance with the criteria of Safe Work Australia - Globally Harmonised System (GHS)

GHS Classification

Reproductive toxicity	Category 1B
-----------------------	-------------

Label elements

Health hazard

**Signal word**

Danger

Hazard statements

May damage the unborn child

Precautionary Statements - Prevention

Obtain special instructions before use

Do not handle until all safety precautions have been read and understood

Use personal protective equipment as required

Precautionary Statements - Response

IF exposed or concerned: Get medical advice/attention

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Other hazards which do not result in classification

May be harmful if inhaled.

Causes mild skin irritation.

Section 3: Composition/information on ingredients

Chemical name	CAS No	Weight-%
Formamide	75-12-7	<70
Non-hazardous ingredients	Proprietary	Balance

Section 4: First aid measures**Description of first aid measures**

General advice	Show this safety data sheet to the doctor in attendance.
Emergency telephone number	Poisons Information Centre, Australia: 13 11 26
Inhalation	Remove to fresh air.
Eye contact	Rinse thoroughly with plenty of water, also under the eyelids.
Skin contact	Wash skin with soap and water.
Ingestion	Rinse mouth.

Most important symptoms and effects, both acute and delayed

Symptoms Prolonged contact may cause redness and irritation.

Indication of any immediate medical attention and special treatment needed

Note to doctors Treat symptomatically.

Section 5: Firefighting measures

Suitable Extinguishing Media

Suitable Extinguishing Media Dry chemical, CO₂, alcohol-resistant foam or water spray.

Unsuitable extinguishing media High volume water jet.

Specific hazards arising from the chemical

Specific hazards arising from the chemical No information available.

Hazardous combustion products Thermal decomposition can lead to release of irritating gases and vapours: Carbon oxides. Nitrogen oxides (NO_x). Sodium oxides. Hydrogen cyanide. Ammonia.

Special protective actions for fire-fighters

Special protective equipment and precautions for fire-fighters Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

Section 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures

Personal precautions Ensure adequate ventilation. Keep people away from and upwind of spill/leak. Evacuate personnel to safe areas.

For emergency responders Use personal protection recommended in Section 8.

Environmental precautions

Environmental precautions Avoid release to the environment. See Section 12 for additional Ecological Information.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so. Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder, sawdust).

Methods for cleaning up Pick up and transfer to properly labelled containers.

Precautions to prevent secondary hazards

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

Section 7: Handling and storage

Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Remove contaminated clothing and shoes. Wear personal protective equipment. Wash hands thoroughly after handling.

General hygiene considerations Do not eat, drink or smoke when using this product. Wash hands before breaks and immediately after handling the product.

Conditions for safe storage, including any incompatibilities

Storage Conditions Store in a well-ventilated place. Store locked up.

Incompatible materials Acids. Bases. Strong oxidising agents. Sulphur trioxide.

Section 8: Exposure controls and personal protection

Control parameters

Exposure Limits

Chemical name	Australia	New Zealand	ACGIH TLV
Formamide 75-12-7	TWA: 10 ppm TWA: 18 mg/m ³	TWA: 10 ppm TWA: 18 mg/m ³ Skin	TWA: 1 ppm S*

Chemical name	European Union	United Kingdom	Germany MAK
Formamide 75-12-7	-	TWA: 20 ppm TWA: 37 mg/m ³ STEL: 30 ppm STEL: 56 mg/m ³	*

Biological occupational exposure limits This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

Appropriate engineering controls

Engineering controls Showers
Eyewash stations
Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).

Skin and body protection Wear suitable protective clothing.

Hand protection Wear suitable gloves: Nitrile rubber.

Respiratory protection No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.

Environmental exposure controls No information available.

Thermal hazards No information available.

Section 9: Physical and chemical properties

Information on basic physical and chemical properties

Appearance

Physical state Liquid
Colour Varies

Odour Odourless
Odour threshold No information available

<u>Values</u>		<u>Remarks • Method</u>
pH		Not applicable
Melting point / freezing point		No data available
Initial boiling point and boiling range		No data available
Flash point	154 °C	
Evaporation rate		No data available
Flammability		Not applicable
Flammability Limit in Air		
Upper flammability or explosive limits		Not applicable
Lower flammability or explosive limits		Not applicable
Vapour pressure		No data available
Vapour density		No data available
Relative density		No data available
Water solubility		No data available
Solubility(ies)		No data available
Partition coefficient		No data available
Autoignition temperature		No data available
Decomposition temperature		No data available
Kinematic viscosity		No data available
Dynamic viscosity		No data available
Explosive properties	No information available.	
Oxidising properties	No information available.	

<u>Other information</u>	
Softening point	No information available
Molecular weight	No information available
VOC Content (%)	No information available
Liquid Density	No information available
Bulk density	No information available
Particle characteristics	No information available

Section 10: Stability and reactivity

Reactivity

Reactivity None under normal use conditions.

Chemical stability

Stability Stable under normal conditions.

Explosion data

Sensitivity to mechanical impact None.

Sensitivity to static discharge None.

Possibility of hazardous reactions

Possibility of hazardous reactions None under normal processing.

Conditions to avoid

Conditions to avoid Excessive heat.

Incompatible materials

Incompatible materials Acids. Bases. Strong oxidising agents. Sulphur trioxide.

Hazardous decomposition products

Hazardous decomposition products Carbon oxides. Nitrogen oxides (NOx). Sodium oxides. Hydrogen cyanide. Ammonia.

Section 11: Toxicological information

Acute toxicity

Information on likely routes of exposure

Product Information

Inhalation May be harmful if inhaled.

Eye contact May cause slight eye irritation.

Skin contact Causes mild skin irritation.

Ingestion Specific test data for the substance or mixture is not available.

Symptoms Prolonged contact may cause redness and irritation.

Numerical measures of toxicity - Product Information

Numerical measures of toxicity No information available.

The following values are calculated based on chapter 3.1 of the GHS document:

ATEmix (oral) 8,022.80 mg/kg

ATEmix (dermal) 8,695.70 mg/kg

ATEmix (inhalation-dust/mist) 30.40 mg/l

Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Formamide	= 5577 mg/kg (Rat)	= 6 g/kg (Rabbit)	> 21 mg/L (Rat) 4 h

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation Classification based on data available for ingredients. May cause skin irritation.

Serious eye damage/eye irritation No information available.

Respiratory or skin sensitisation No information available.

Germ cell mutagenicity No information available.

Carcinogenicity No information available.

Reproductive toxicity Contains a known or suspected reproductive toxin. Classification based on data available for ingredients. May damage fertility or the unborn child.

STOT - single exposure No information available.
STOT - repeated exposure No information available.
Aspiration hazard No information available.

Section 12: Ecological information

Ecotoxicity

Aquatic ecotoxicity The environmental impact of this product has not been fully investigated.

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Formamide	EC50: >500mg/L (72h, Desmodesmus subspicatus) EC50: >500mg/L (96h, Desmodesmus subspicatus)	LC50: =9135mg/L (96h, Brachydanio rerio)	-	EC50: >500mg/L (48h, Daphnia magna)

Terrestrial ecotoxicity There is no data for this product.

Persistence and degradability

Persistence and degradability No information available.

Bioaccumulative potential

Bioaccumulation

Component Information

Chemical name	Partition coefficient
Formamide	-0.82

Mobility

Mobility No information available.

Other adverse effects

Other adverse effects No information available.

Section 13: Disposal considerations

Waste treatment methods

Waste from residues/unused products Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

Contaminated packaging Do not reuse empty containers.

See section 8 for more information

Section 14: Transport information

ADG Not regulated
IATA Not regulated
IMDG Not regulated

Transport in bulk according to Annex II of MARPOL and the IBC Code
 No information available

Section 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulations

Australia

Classified as a hazardous substance in accordance with the criteria of Safe Work Australia - Globally Harmonised System (GHS).

See section 8 for national exposure control parameters

Standard for the Uniform Scheduling of Medicines and Poisons (SUSMP)

No poisons schedule number allocated

Poison Schedule Number Not applicable

Australian Industrial Chemicals Introduction Scheme (AICIS)

Contact supplier for inventory compliance status

Illicit Drug Precursors/Reagents

This product does not contain any substance(s) on the Illicit Drug Precursors/Reagents list.

Chemical name	Illicit Drug Precursors/Reagents
Formamide - 75-12-7	Category 2

National pollutant inventory

Subject to reporting requirement

Chemical name	National pollutant inventory
Formamide - 75-12-7	20 MW Threshold category 2b total 60000 MWH Threshold category 2b total 1 tonne/h Threshold category 2a total 25 tonne/yr Threshold category 1a total 400 tonne/yr Threshold category 2a total 2000 tonne/yr Threshold category 2b total

International Inventories

AICS Contact supplier for inventory compliance status.
NZIoC Contact supplier for inventory compliance status.
TSCA Contact supplier for inventory compliance status.
DSL/NDSL Contact supplier for inventory compliance status.

EINECS/ELINCS	Contact supplier for inventory compliance status.
ENCS	Contact supplier for inventory compliance status.
IECSC	Contact supplier for inventory compliance status.
KECL	Contact supplier for inventory compliance status.
PICCS	Contact supplier for inventory compliance status.

Legend:

AICS	- Australian Inventory of Chemical Substances
NZIoC	- New Zealand Inventory of Chemicals
TSCA	- United States Toxic Substances Control Act Section 8(b) Inventory
DSL/NDL	- Canadian Domestic Substances List/Non-Domestic Substances List
EINECS/ELINCS	- European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
ENCS	- Japan Existing and New Chemical Substances
IECSC	- China Inventory of Existing Chemical Substances
KECL	- Korean Existing and Evaluated Chemical Substances
PICCS	- Philippines Inventory of Chemicals and Chemical Substances

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

Section 16: Any other relevant information

Issuing Date	14-Oct-22
Revision Date	14-Oct-22
Revision Note	Initial Release.

Key or legend to abbreviations and acronyms used in the safety data sheet**Legend Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION**

TWA	TWA (time-weighted average)	STEL	STEL (Short Term Exposure Limit)
Ceiling	Maximum limit value	*	Skin designation
C	Carcinogen		

Key literature references and sources for data used to compile the SDS

Agency for Toxic Substances and Disease Registry (ATSDR)
 U.S. Environmental Protection Agency ChemView Database
 European Food Safety Authority (EFSA)
 EPA (Environmental Protection Agency)
 Acute Exposure Guideline Level(s) (AEGl(s))
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
 U.S. Environmental Protection Agency High Production Volume Chemicals
 Food Research Journal
 Hazardous Substance Database
 International Uniform Chemical Information Database (IUCLID)
 Japan GHS Classification
 Australian National Industrial Chemicals Notification and Assessment Scheme (NICNAS)
 Australian Industrial Chemicals Introduction Scheme (AICIS)
 NIOSH (National Institute for Occupational Safety and Health)
 National Library of Medicine's ChemID Plus (NLM CIP)
 National Library of Medicine's PubMed database (NLM PUBMED)
 National Toxicology Program (NTP)
 New Zealand's Chemical Classification and Information Database (CCID)

Organisation for Economic Co-operation and Development Environment, Health, and Safety Publications
Organisation for Economic Co-operation and Development High Production Volume Chemicals Programme
Organisation for Economic Co-operation and Development Screening Information Data Set
World Health Organization

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

End of Safety Data Sheet